1. Uso Previsto

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14. Bibliografia

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7): 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1): 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3): 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5): 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

Rev. 1, 2014-09-12 AF

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limitación de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour pour "n" tests	Válido para para <n> ensayos	Sufficiente per "n" saggi	Wystarczający na <n> powtórzeń

Trade diagnostical products
for clinical laboratories
products for life.

Italy Corporate Headquarter
Meridian Healthcare
Via G. Guglielmino, 65
95030 Trimestri Estero (Catania) – Italy
Tel: +39 095 7256869 / +39 095 419 88 60
Fax: +39 095 725 44

Customer Service
Meridian Healthcare
Tel: +39 095 419 88 60 / +39 0922 83 21 88
Fax: +39 095 725 44 54
mail: info@meridianhealthcare.it
marketing@meridianhealthcare.it

Technical Support
Engineering Department
Meridian Healthcare
Tel: +39 095 419 88 61
Fax: +39 095 725 44 54
mail: specialist@meridianhealthcare.it

Sales and Orders
Meridian Healthcare
Tel: +39 095 725 68 69
Fax: +39 095 725 44 54
mail: ordini@meridianhealthcare.it
amministrazione@meridianhealthcare.it