

Quantitative determination of Phenytoin
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

This is a quantitative turbidimetric test for the measurement of Phenytoin in human serum.

Latex particles coated with Phenytoin are agglutinated when mixed with Phenytoin Antibody solution. When a sample containing Phenytoin is used, the agglutination reaction is partially inhibited, slowing the process of agglutination. The agglutination causes an absorbance change, inversely dependent upon the Phenytoin contents of the patient sample.

CLINICAL SIGNIFICANCE

Phenytoin is an antiepileptic drug that is used to control focal and generalized motor seizures (fits, tonic-clonic). The range of serum concentrations at which Phenytoin toxicity occurs is broad. Toxic effects, such as nystagmus and distaxia can occur at concentrations above 20 µg/mL. Ataxia and eye movement abnormalities may occur above 30 µg/mL.

REAGENTS

R 1 Antibody Buffer	Anti- Phenytoin antibody. Sodium azide 0.09% w/v
R2 Latex Reagent	Latex. Sodium azide 0.09% w/v

CALIBRATION

The use of Therapeutic Drug Calibrator is recommended. We recommend a Multi Point calibration each 7 days, each change of reagent bottle/ lot or according to quality control procedures.

PREPARATION

Reagents: Ready for use.

Calibrator/Control: Lyophilized. Reconstitute with the indicated volume on the respective instructions of use.

Calibration Curve: Must be used the Spinreact Therapeutic Drug Calibrator.

STORAGE AND STABILITY

R1. Antibody Buffer

Stable up to expiry date when stored at 2-8°C. Do not freeze.

Before use, reagents should be gently swirled to dislodge bubbles and ensure homogeneity.

R2. Latex Reagent

Stable up to expiry date when stored at 2-8°C. Do not freeze.

Before use, reagents should be gently swirled to dislodge bubbles and ensure homogeneity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 700 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Samples should be stored at between 2-8°C prior to testing for up to three days. Samples intended for assaying after three days from collection should be frozen at -20°C until use. Any additional clotting or precipitation that occurs due to the freeze/thaw treatment should be removed by centrifugation prior to analysing the Phenytoin concentration of that sample.

PROCEDURE

1. Bring the reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions:
 - Wavelength : 700
 - Temperature : 37 °C
 - Cuvette ligh path : 1cm
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette:

Reagent R1	990 µL
Sample or Calibrator	9 µL
Reagent R2	360 µL

5. Mix and read the absorbance (A1) after the R2 addition.

6. Incubate at 37°C and read the absorbance (A2) exactly 4 minutes after the R2 addition.

CALCULATIONS

Calculate the absorbance difference (A2-A1) of each point of the calibration curve and plot the values obtained against the Phenytoin concentration of each calibrator dilution. Phenytoin concentration in the sample is calculated by interpolation of its (A2-A1) in the calibration curve.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact has the Therapeutic Drug Controls kit Ref.:939550 (3 levels). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

10 – 20 µg/mL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Assay Range: The range of this assay is approximately 1.1- 40 µg/mL depending on the concentration range of the Phenytoin calibrators in use.

Samples with concentrations in excess of the highest calibrator should be diluted with the 0 µg/mL calibrator, re-assay and multiply the result with the appropriate dilution factor.

Precision:

	Intra-assay (n=20)			Inter-assay (n=20)		
Mean (µg/mL)	4.99	15.80	25.01	4.52	14.55	21.30
SD	0.184	0.268	0.461	0.220	0.340	0.468
CV (%)	3.68	1.69	1.84	4.88	2.34	2.20

Accuracy: Results obtained using reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 40 samples spanning the range 2.37 to 39.3 µg/mL were the following:

Correlation coefficient (r): 0.9962.

Regression equation: $y = 1.0515x + 1.635$

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

The following analytes were tested up to the following levels and were found not to interfere:

Intralipid®	800 mg/dl
Bilirubin	25 mg/dl
Triglyceride	1000 mg/dl
Haemoglobin	1000 mg/dl

NOTES

1. In order to avoid contamination it is recommended to use disposable material.
2. **M.H. has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Kodama Y., Koike Y., Kimoto H., Yasuraga F., Takeyama M., Teraoka I., Fujii I.: Binding parameters of valproic acid to serum protein in healthy adults at steady state. *Therapeutic Drug Monitoring* 1992, **14**; 55-60.
2. Tisdale J.E., Tsuyuki R. T., Oles K. S., Penry J.K.: Relationship between serum concentrations and dose of valproic acid during monotherapy in adult outpatients. *Therapeutic Drug Monitoring* 1992, **14**; 416-423.
3. Chadwick D.W.: Valproate monotherapy in the Management of generalized and partial seizures. *Epilepsia* 1987, **78 (supplement 2)**: 512-7.
4. Chadwick D.W.: Concentration effect relationships of Valproic Acid. *Clinical Pharmacokinetics* 1985, **10**; 155-163.
5. Newman, D.J., Henneberry H., Price C.P., "Particle Enhanced Light Scattering Immunoassay", *Ann. Clin. Biochem.*, **29**: 22-42 (1992)

REF. 939050

Cont.

2 x 17 mL, 2 x 6 mL


Meridian Healthcare srl

 Tel. +39 095 725 68 69 Fax: +39 095 725 44 54
 info@meridianhealthcare.it
 www.meridianhealthcare.it
