
Capilia™ hMPV
Metapneumovirus Ag
Rapid Test

CAPILIA™ METAPNEUMOVIRUS
LATERAL FLOW RAPID TEST

To detect human metapneumovirus antigens in nasal swab, nasal aspirate or pharyngeal swab (to assist in the diagnosis of human metapneumovirus infection)

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Capilia™ METAPNEUMOVIRUS

Rapid Test

TO DETECT HUMAN METAPNEUMOVIRUS ANTIGENS IN NASAL SWAB, NASAL ASPIRATE OR PHARYNGEAL SWAB (TO ASSIST IN THE DIAGNOSIS OF HUMAN METAPNEUMOVIRUS INFECTION)

Human metapneumovirus (hMPV), first identified in the Netherlands in 2001, is an RNA virus that belongs to the *Metapneumovirus* genus within the *Pneumovirinae* subfamily of the *Paramyxoviridae* family. hMPV is a virus that can cause respiratory infection. hMPV infection is a typical respiratory infection in infants. It is similar to respiratory syncytial (RS) virus infection and differentiation from each other is difficult. hMPV infection is an important factor for inducing bronchiolitis. The majority experience primary hMPV infection during infancy, and re-infection occurs frequently. As hMPV infection may become prevalent in infants and the elderly and increase in severity, early detection is important.

Conventionally, genetic testing is used for the detection of hMPV, but it has some disadvantages, for instance, requiring complicated manipulation, special equipment and instruments, and a certain period of time to obtain the test results. Compared to genetic testing, this product allows rapid detection of hMPV, without requiring special skills or instruments.

PRINCIPLE OF THE TEST

Measurement using this product is based on an immunochromatography assay using a monoclonal antibody that recognizes human metapneumovirus antigens.

This product comprises a test plate with a carrier strip containing a sample placement area, a reagent area including a colloidal platinum-gold labeled anti-human metapneumovirus monoclonal antibody (mouse) (hereinafter referred to as "colloidal platinum-gold labeled anti-human metapneumovirus antibody"), and a developing area that fixes an anti-human metapneumovirus monoclonal antibody (mouse) (hereinafter referred to as "anti-human metapneumovirus antibody") and an anti-mouse immunoglobulin polyclonal antibody (rabbit) (hereinafter referred to as "anti-mouse immunoglobulin antibody").

When a sample is placed on the sample placement area of the test plate, the colloidal platinum-gold labeled anti-human metapneumovirus antibody dissolves and forms an immune complex with human metapneumovirus antigens in the sample. This immune complex migrates through the developing area by capillary action, is captured by the anti-human metapneumovirus antibody fixed in the developing area, and forms a black line of colloidal platinum-gold in the reading area [T]. The black line visually displays the existence of human metapneumovirus antigens in the sample.

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Nasal Swab

Correlation with PCR method:

- sensitivity: 92,0%
- specificity: 100%
- accuracy: 96,3%

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
Genotipo B: $3,0 \times 10^2$ TCID₅₀/test

DETECTION LIMIT

The minimum detection limit of this product is as follows:

Genotype A: $4,5 \times 10^2$ TCID₅₀/test

Genotype B: $3,0 \times 10^2$ TCID₅₀/test

INFO ORDINI / ORDERING INFORMATION

Capilia™ hMPV

REF: CAHM1670

N. TEST: 20