

ACL® ANALYZER ROTORS

Cuvette Rotor for use on IL ACL® Analyzer models: 3000, 7000, 9000 & ELITE-PRO



- Excellent plastic and absolute transparency for a reading of the sample to 660 and 405 nm.
- 100% overlap with the original rotors
 - 100 Rotors
 - 20 Position
- Controlled Laser system one by one.

ORDERING INFORMATION

REF.	PRODUCT	PACK-SIZE
ACL-AR-I-100	ACL ANALYZER ROTORS	100 pcs

ACL® ANALYZER ROTORS

Rotori per Analizzatori ACL®, serie 1000, 2000, 3000, 6000, 7000, 9000 & ELITE-PRO

- Plastica altamente trasparente per una lettura ottimale sui canali coagulativi e cromogenici a 660 e 405 nm.
 - 100% sovrapponibili agli originali
 - 100 Rotori
 - 20 Posizioni
 - Controllati al Laser singolarmente

