

ADENOVIRUS PNEUMONIAE Rapid Test

The Adenovirus Antigen Rapid Test is for the qualitative detection of Adenovirus antigen in eye conjunctive swab, throat swab and nasal swab as an aid in the diagnosis of adenovirus infections.

FEATURES

- Simple, Accurate, Fast.
- Eye swab, throat swab and nasal swab.
- 🕒 Reliable Results in 15 minutes.
- Sensitivity: 98.6%
- Specificity: 98.2%
- Accuracy: 95.6%
- In Vitro Diagnostics



ORDERING INFORMATION:

Code	Product	Pack-size
IAD-502	ADENO Pneum. Rapid test	20 test



ADENOVIRUS PNEUMONIAE Rapid Test

Test immunocromatografico Lateral Flow per la rilevazione dell'ADENOVIRUS respiratorio in secrezioni naso-faringee e da tampone oculare.

SPECIFICHE

- Semplice, Accurato, Veloce.
- Tampone Nasofaringeo e Oculare.
- 🕒 Risultato in 15 minuti.
- Sensibilità: 98.6%
- Specificità: 98.1%
- Accuratezza: 98.2%
- Per Uso Diagnostico in Vitro

