

PANKREAS α-AMYLASE EPS FLUID IFCC

Enzymatic colorimetric determination of pancreatic alpha-amylase according to the IFCC in serum and plasma. For *in vitro* diagnostic use only.



- Enzymatic Colorimetric Assay.
- Dual Vial Liquid Stable Test Kit.
- ⌚ Reliable Results in 8 minutes.
- Measuring range: 2 - 1300 U/L.
- Serum or Plasma acceptable.
- Applications are available for many analyzers after request.



ORDERING INFORMATION:



Product Code
AF06000030

Product Open System
Amylase-EPS-Pancr. IFCC fluid

Pack-Size
2x25 + 2x5,5 ml.

α-AMILASI PNCREATICA EPS FLUID IFCC

Saggio enzimatico colorimetrico per la determinazione dell'Amilasi Pancreatica su campioni di siero, plasma e urine.

- Metodo Enzimatico Colorimetrico IFCC.
- Due reagenti liquidi.
- ⌚ Risultato in 8 minuti.
- Saggio su siero, plasma e urine.
- Automatizzabile su tutte le piattaforme open system.
- Range di misura: 2 - 1300 U/L.

DEDICATED PACKAGING

Product Code

AF06911060
AF06911120-B
AF06917084-B
AF06917084
AF06SELO30

ANALYZER

HITACHI 912
HITACHI 911 BARCODED
HITACHI 917 BARCODED
OLYMPUS AU
SELECTRA SERIES

Pack-Size

6x50 + 3x20 ml.
3x100 + 3x20 ml.
3x67 + 3x14 ml.
3x67 + 3x14 ml.
5x25 + 5x5 ml.

