

BIOPANDA™ CELIAC DISEASE tTg IgA Rapid Test

The **Celiac Disease tTG** rapid test is an immunochromatographic test designed to detect, in human blood, IgA-type antibodies against human tissue transglutaminase (tTG), the main autoantigen recognized by the anti-endomysial antibodies.



- Direct Testing in serum, plasma and wb.
 - Sample: 10µL. s/p and 20µL. wb
 - ⌚ Reliable Results in 10 minutes.
 - Sensitivity tTG : 95,0%
 - Specificity tTG: 98,5%
 - Long shelf life of 28 months



BIOPANDA™ CELIAC DISEASE tTg IgA Test Rapido

Il dispositivo **Celiachia BIOPANDA™ tTG** è un test rapido Lateral Flow per la determinazione degli anticorpi anti-Tranglutaminasi IgA su campioni di sangue intero, siero o plasma. L'utilizzo è particolarmente indicato nel caso di screening in laboratorio.

- Semplice, Accurato, Veloce.
- Test diretto su sangue intero, siero o plasma
- Campione: 10µL. per s/p e 20µL per sangue venoso o capillare.
- ⌚ Risultato in 10 minuti.
- Sensibilità tTG : 95,0%
- Specificità tTG: 98,5%



ORDERING INFORMATION

Codice Dispositivo
RAPG-CEL-001

Dispositivo
BIOPANDA™ Celiac tTG IgA Rapid Test

Confezione
20 test