

# CALPROTECTIN Rapid Test

RAPID-VIDITEST Calprotectin is a rapid test intended for the qualitative detection of calprotectin in human feces specimens, which might be useful for the diagnosis of inflammatory gastrointestinal disorders. This assay is useful in differentiating organic (IBD) from functional gastrointestinal disease (Intestinal Bowel Syndrome). It is a simple, non-invasive biomarker that is especially useful in children, who may require general anesthesia for colonoscopy. Only for laboratory and gastroenterologist use.

- ⌚ Simple, Accurate, Fast.
- Direct Testing on human feces specimens.
- Ambulatory can be used.
- ⌚ Reliable Results in 10 minutes.
- Sensitivity : > 94,0%
- Specificity: > 93,0%
- Cut-off: 50 µg hCp/g feces.
- CE marked



# CALPROTECTINA Rapid Test

Test rapido immunocromatografico lateral flow per la rilevazione qualitativa della Calprotectina in campioni fecali e per differenziare l'IBD dalla malattia intestinale funzionale (IBS). La membrana è pretrattata nella zona test con anticorpi monoclonali anti-calprotectina.



- ⌚ Semplice, Accurato, Veloce.
- Test diretto su campioni di feci.
- Utilizzo previsto anche in ambulatorio.
- ⌚ Risultato in 10 minuti.
- Sensibilità: > 94,0%.
- Specificità: > 93,0%
- Cut-off: 50 µg hCp/g feci.
- CE marked

EUROPE Ordering Information:  
CALPROTECTIN RAPID VIDITEST  
Product Code: ODZ-248  
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 -725 68 69 (0039)