

CAMPYLOBACTER SPECIES Ag Rapid Test

The Rapid-VIDITEST Campylobacter is a rapid chromatographic immunoassay for the qualitative detection of Campylobacter in feces specimens, which might be useful for the diagnosis of campylobacteriosis. Only for laboratory use.



- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- 🕒 Reliable Results in 10 minutes
 - Sensitivity : >99%
 - Specificity: >99%
- Campylobacter jejuni
- Campylobacter coli
 - CE marked

CAMPYLOBACTER SPECIES Ag Rapid Test

Test rapido per la determinazione qualitativa di Campylobacter spp.(identifica le specie patogene Campylobacter jejuni e Campylobacter coli) in campioni fecali, utile nella diagnosi di campylobatteriosi

- Semplice, Accurato, Veloce
- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 10 minuti
- Sensibilità : >99%
- Specificità: >99%
- Campylobacter jejuni
- Campylobacter coli
- CE marked

Controlli disponibili:

- Ctrl Positivo: ODZ-208 kit da 10 test
- Ctrl Positivo: RM01340PC fiala da 1 ml



EUROPE Ordering Information:
CAMPYLOBACTER RAPID TEST

Product Code: ODZ-160

Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)