


# CANDIDA Albicans YEAST



Vaginal Swabs & URINE


## CANDIDA Albicans YEAST Rapid Test

New Rapid Test Lateral Flow for the qualitative detection of CANDIDA Albicans in swab from vaginal cavity and URINA.

- Simple, Accurate, Fast
- Sample: vaginal swab and URINA
- ⌚ Reliable Results in 15 minutes
- Sensitivity: 95,5%
- Specificity: 98,4%

### PERFORMANCE

The gold standard was the patient positive to speculum and vaginal examinations to be diagnosed Candida Albicans Yeast infection. The total 300 cases were compared to wet mount microscopy pouch culture and JD's Candida Kit.  
**Predictive Value JD's Candida Albicans YEAST:**  
**Positive 97,2%**  
**Negative 97,4%**



### ORDERING INFORMATION:

**Product Code**  
JCAAG-02

**Product**  
CANDIDA Albicans Rapid Test

**Pack Size**  
25 test

