

CRYPTOSPORIDIUM PARVUM Rapid Test

The Cryptosporidium VIDITEST® is a rapid chromatographic immunoassay for the qualitative detection of *Cryptosporidium* antigens in human feces specimens to aid in the diagnosis of cryptosporidiosis. Only for laboratory use.



Features

- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- ⌚ Reliable Results in 10 minutes
 - Sensitivity : >99%
 - Specificity: >99%
 - CE marked



CRYPTOSPORIDIUM PARVUM Test Rapido

Test immunocromatografico in vitro per la determinazione degli oocisti di *Cryptosporidium parvum* in campioni di feci.

Specifiche

Semplice, Accurato, Veloce

- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- ⌚ Risultato in 10 minuti
- Sensibilità : >99%
- Specificità: >99%
- Limite di risoluzione: 50-100 oocisti / 100µl (diluizione di un campione positivo mediante metodo di Ziels-Niehlsen).
- CE marked



Controlli disponibili:

- Ctrl Positivo: ODZ-205 kit da 10 test



EUROPE Ordering Information:
CRYPTOSPORIDIUM PARVUM RAPID TEST
Product Code: ODZ-278
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)

