

ESCHERICHIA COLI O157:H7 Rapid Test

The Rapid-VIDITEST *E.coli* O157:H7 Card is a one step coloured chromatographic immunoassay for the qualitative detection of *Escherichia coli* O157 in stool samples. Only for laboratory use.



Features

- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- 🕒 Reliable Results in 10 minutes
 - Sensitivity : 100%
 - Specificity: 85%
 - CE marked



ESCHERICHIA COLI O157:H7 Test Rapido

O157 E. COLI CARD è un test immunocromatografico manuale e rapido per la rilevazione qualitativa di *Escherichia coli* O157:H7 negli alimenti e nei campioni di feci umane, di ausilio nella diagnosi delle infezioni da *E. coli*. È possibile rilevare l'*E. coli* O157:H7 da colonie sospette in colture di feci. Il test offre un'analisi di screening semplice e altamente sensibile per effettuare una diagnosi presuntiva di infezione da *Escherichia coli* O157:H7 e può essere utilizzato per identificare isolati sospetti di *E. coli* O157:H7 da terreni selettivi.

Specifiche

- Semplice, Accurato, Veloce
- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 10 minuti
- Sensibilità : 100%
- Specificità: 85%
- CE marked

Controlli disponibili:

- Ctrl Positivo: ODZ-207 kit da 10 test



EUROPE Ordering Information:
ESCHERICHIA COLI O157:H7 RAPID TEST
Product Code: ODZ-187
Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)