

## FAECAL OCCULT BLOOD (FOB) Rapid Test

The M.H. Faecal Occult Blood (FOB) Test Device is a rapid, visual immunochromatographic test for the qualitative detection of human haemoglobin in faecal samples. The test is for professional and in vitro use only. This test was created to detect indications of gastrointestinal (g.i.) disorders like bowel cancer and major haemorrhoids.

- 🕒 Simple, Accurate, Fast.
- Direct Testing on human feces specimens.
- Ambulatory can be used.
- 🕒 Reliable Results in 5 minutes.
- Sensitivity : 100%
- Specificity: 99,9%
- Cut-off: 40 ng/ml.
- CE marked



## SANGUE OCCULTO nelle FECEI Test Rapido

Hb Fecale è un test rapido, immunocromatografico visivo per la determinazione qualitativa dell'emoglobina umana in campioni di feci. Questo test vuole essere un aiuto nella diagnosi di malattie del basso tratto gastrointestinale. Il test è raccomandato per uso professionale.



- 🕒 Semplice, Accurato, Veloce.
- Test diretto su campioni di feci.
- Utilizzo previsto anche in ambulatorio.
- 🕒 Risultato in 5 minuti.
- Sensibilità: 100%.
- Specificità: 99,9%
- Cut-off: 40 ng/ml.
- CE marked

### EUROPE Ordering Information:

**FOB RAPID VIDITEST**

**Product Code: 272001**

**Packaging: 20 test**



Meridian Healthcare®

Meridian Healthcare srl

68 - Tremestieri Etneo

Catania, IT 95030

Phone: 095 -725 68 69 (0039)