

FOOD INTOLERANCE

90 Foods IgG ELISA Kit

For Semi-Quantitative Analysis of Antibodies to 90 Food Allergens in Human Serum

SPECIFICHE

- ❖ Mediterranean Panel
- ❖ Food-Specific IgG antibody
- ❖ 3 plates x 96 wells/plate ELISA Kit
- ❖ Foods IgG Positive Control included
- ❖ Sensitivity : 99%
- ❖ Specificity: 99%
- ❖ Reproducibility: 100%



90 FOODS SPLIT MICROPLATE MAP

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK	Apple	Broccoli	Chard	Codfish	Garlic	Lattuce, Iceberg	Mustard Seed	Pear	Rice	Squashes	Trout
B	Calibrator 1	Artichoke	Butter	Cheese (cottage)	Coffee	Grape, White/con cord	Lemon	Oat	Black Pepper	Rye	Squid	Tuna
C	Calibrator 2	Asparagus	Cabbage	Cheese (Cured)	Cola nut	Grapefruit	Lentils	Olive	Pinto bean	Salmon	Strawberry	Turkey
D	Calibrator 3	Avocado	Cane Sugar	Chick peas	Corn	Green pea	Lima bean	Onion	Pineapple	Sardine	String bean	Walnut, black
E	Calibrator 4	Banana	Cantaloupe	Chicken	Cow's milk	Green pepper	Lobster	Orange	Plum	Shrimp	Sunflower seed	Wheat
F	Positive Control	Barley, whole grain	Carrot	Chocolate	Cucumber	Hake	Malt	Parsley	Pork	Sole	Sweet Potato	Yeast, Baker's
G	Almond	Beef	Cauliflower	Cinnamon	Egg, white/yolk	Honey	Marjoram	Peach	Potato	Soybean	Tea, Black	Yeast, Brewer's
H	American Cheese	Beets	Celery	Clam	Eggplant	Lamb	Mushroom	Peanut	Rabbit	Spinach	Tomato	Yogurt

Meridian Healthcare srl

Tel. +39 095 725 68 69 Fax: +39 095 725 44 54

info@meridianhealthcare.it

www.meridianhealthcare.it

Meridian Healthcare®
medical laboratory equipment

Meridian Healthcare®

