



GARDNERELLA Vaginalis



vaginal swab
glans
URINA




GARDNERELLA Vaginalis Rapid Test

New Rapid Test Lateral Flow for the qualitative detection of Gardnerella Vaginalis in swab from vaginal cavity, Glans and URINA.

- Simple, Accurate, Fast
- Sample: vaginal swab, glans and URINA
- ⊕ Reliable Results in 15 minutes
- Sensitivity: 98,46%
- Specificity: 98,62

PERFORMANCE

Due to the nature of the test system, precision was evaluated by testing various samples on consecutive days, by multiple clinicians at different geographic locations. JD's GardnerellaV. antigen Test precision studies were run on artificially generated laboratory samples reflecting the known range of Gardnerellavaginalis parasite burden levels in female patients. Given the design of the product, the JD's GardnerellaV. antigen Test performed according to expectations: the JD's GardnerellasV. antigen test kit performed with >99% precision over the intra-assay (between day, betweenclinician, and between clinical site) parameters.



ORDERING INFORMATION:

Product Code
JGVAG-02

Product
GARDNERELLA Vaginalis Rapid Test

Pack Size
25 test

