

GIARDIA LAMBLIA Ag Rapid Test

The Rapid-VIDITEST Giardia is a rapid chromatographic LATERAL FLOW immunoassay for the qualitative detection of Giardia antigens in human feces specimens to aid in the diagnosis of giardiasis. Only for laboratory use.



- Simple, Accurate, Fast
- Direct Testing on Faeces
- Ambulatory can be used
- 🕒 Reliable Results in 10 minutes
 - Sensitivity : >97%
 - Specificity: >99%
 - CE marked

No-cross reactivity with the following intestinal pathogens:
Cryptosporidium parvum, Entamoeba histolytica

GIARDIA LAMBLIA Ag Rapid Test

Test rapido immunocromatografico Lateral Flow per la determinazione qualitativa dell'antigene Giardia Lamblia in campioni fecali.

- Semplice, Accurato, Veloce
- Test diretto su campioni fecali
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 10 minuti
- Sensibilità : >97%
- Specificità: >99%
- CE marked

Nessuna cross-reattività con altri patogeni intestinali, tipo:
Cryptosporidium parvum, Entamoeba histolytica

Controlli disponibili:

- Ctrl Positivo: ODZ-206 kit da 10 test



EUROPE Ordering Information:

GIARDIA LAMBLIA RAPID TEST

Product Code: ODZ-016

Packaging: 20 test



Meridian Healthcare®

Meridian Healthcare srl

68 - Tremestieri Etneo

Catania, IT 95030

Phone: 095 725 6869 (0039)

Fax: 095 725 4454 (0039)