

# HELICOBACTER Pylori Ab Rapid Test

JusChek® H. pylori Antibody Rapid Test is a rapid visual immunoassay for the qualitative presumptive detection of specific IgM and IgG antibodies to Helicobacter Pylori with human whole blood, serum and plasma as specimen.




- Direct Testing on serum, plasma and whole blood.
  - Fingertip blood can be used.
- ⌚ Reliable Results in 10 minutes.
  - Sensitivity : 96,8%
  - Specificity: 93,0%
  - Accuracy: 94,6%
  - CE marked

# HELICOBACTER Pylori Ab Rapid Test

JusChek® H.Pylori test rapido immunocromatografico lateral flow per rilevazione qualitativa degli anticorpi IgG e IgM dell'Helicobacter Pylori su campioni di siero, plasma e sangue intero.

- ⌚ Semplice, Accurato, Veloce.
- Test diretto su campioni di siero, plasma e sangue intero.
- Utilizzo previsto anche in ambulatorio e in pediatria.
- ⌚ Risultato in 10 minuti.
- Sensibilità: 96,8%.
- Specificità: 93,0%
- Accuratezza: 94,6%
- CE marked



EUROPE Ordering Information:  
**HELICOBACTER PYLORI Ab**  
Product Code: IHP-402  
Packaging: 25 test



Meridian Healthcare®

Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 -725 68 69 (0039)