

HEPATITIS DISEASE IMMUNOASSAY LATERAL FLOW



AssureTech® EPATITE A

Nuovo Test Rapido Lateral Flow per lo screening rapido qualitativo degli anticorpi del Virus dell'Epatite A in campioni di siero o plasma. Test professionale per esclusivo uso diagnostico in vitro.

HAV IgM

SPECIFICHE:

- Test diretto su siero o plasma.
- Rilevazione IgM
- 🕒 Tempo di esecuzione: 15 min.
- Campione: 5 µL.
- Confezione: 25 test
- Sensibilità relativa: > 99,9%
- Specificità relativa: 100%
- Accuratezza analitica: 99,6%



AssureTech® HEPATITIS A

HAV IgM Rapid Test is an in vitro diagnostic kit for the qualitative detection of anti-HAV IgM in human, serum or plasma specimens. HAV IgM Rapid Test is intended to be used as an aid in the diagnosis of hepatitis A virus infection.

HAV IgM

PERFORMANCE CHARACTERISTICS:

- Direct Testing from serum or plasma.
- Detects IgM
- 🕒 Reliable Results in 15 minutes.
- Sample: 5 µL.
- Packaging: 25 test
- Relative Sensitivity: > 99,9%
- Relative Specificity: 100%
- Accuracy: 99,6%

ORDERING INFO:

REF. HAV-W23

Product: HAV IgM

Packaging: 25 test

