

Attachment 1 - Allegato 1

(en) Update of PT-HRTF expected values (seconds, %, INR) for Sclavo PT-HRTF Lot 55
 (IT) Aggiornamento valori attesi PT-HRTF (secondi, %, INR) per Sclavo PT-HRTF Lotto 55

(en) SCLAVO NORMAL CONTROL PLASMA KIT
 (it) SCLAVO PLASMA DI CONTROLLO NORMALE

REF K29005301



LOT

55

Kit
 Vial / Flacone



2025/02/22

LEGENDA					
V.T.	VALORE TEORICO			R.A.	RANGE ATTESO
T.V.	TARGET VALUE			E.R.	EXPECTED RANGE

TAVOLA DEI VALORI ATTESI - TABLE OF EXPECTED VALUES

REAGENTE SCLAVO SCLAVO REAGENT	PARAMETRO PARAMETER	IL ACL7000/9000/ELITE PRO			SYSMEX CA-5xx/CA-6xx			SYSMEX CA-1500/CA-7000			SYSMEX CS-1600			SYSMEX CS-2100/CS-2500			SYSMEX CS-5100		
		R.A.	E.R.	V.T.	R.A.	E.R.	V.T.	R.A.	E.R.	V.T.	R.A.	E.R.	V.T.	R.A.	E.R.	V.T.	R.A.	E.R.	V.T.
PT-ReadyQuick	Seconds	10.8	16.2	13.5	9.8	14.8	12.3	9.8	14.8	12.3	9.8	14.8	12.3	9.8	14.8	12.3	9.8	14.8	12.3
	%	87	129	108	89	125	107	89	125	107	89	125	107	89	125	107	89	125	107
	INR	0.78	1.16	0.98	0.82	1.12	0.97	0.82	1.12	0.97	0.82	1.12	0.97	0.82	1.12	0.97	0.82	1.12	0.97
PT-HRTF	Seconds	9.2	13.8	11.5	9.0	13.5	11.2				9.0	13.5	11.2	9.0	13.5	11.2	9.0	13.5	11.2
	%	70	106	88	73	101	87				73	101	87	73	101	87	73	101	87
	INR	0.85	1.27	1.06	0.98	1.28	1.13				0.98	1.28	1.13	0.98	1.28	1.13	0.98	1.28	1.13
aPTT-S	Seconds	25.3	37.9	31.6	21.1	31.1	26.1	21.1	31.1	26.1	21.1	31.1	26.1	21.1	31.1	26.1	21.1	31.1	26.1
Fibrinogen	mg/dL	200	320	260	227	307	267	227	307	267	227	307	267	227	307	267	227	307	267
Antithrombin (FXa)	%	67	101	84	76	96	86	76	96	86	76	96	86	76	96	86	76	96	86
TT	Seconds	14	22	18	14.9	22.9	18.9	17.9	25.9	21.9	12.8	20.8	16.8	12.8	20.8	16.8	12.8	20.8	16.8
Protein C - Chromogenic	%	70	104	87	84	104	94	84	104	94	84	104	94	84	104	94	84	104	94
Protein S - Clotting	%	77	115	96	71	95	83	71	95	83	71	95	83	71	95	83	71	95	83
Free Protein S	%	99	149	124	90	122	106	90	122	106	90	122	106	90	122	106	90	122	106
Factor II	%				75	95	85	75	95	85	75	95	85	75	95	85	75	95	85
Factor V	%				78	98	88	78	98	88	78	98	88	78	98	88	78	98	88
Factor VII	%				83	103	93	83	103	93	83	103	93	83	103	93	83	103	93
Factor VIII	%				98	118	108	98	118	108	98	118	108	98	118	108	98	118	108
Factor IX	%				90	110	100	90	110	100	90	110	100	90	110	100	90	110	100
Factor X	%				77	97	87	77	97	87	77	97	87	77	97	87	77	97	87
Factor XI	%				81	109	95	81	109	95	81	109	95	81	109	95	81	109	95
Factor XII	%				87	117	102	87	117	102	87	117	102	87	117	102	87	117	102



CA-6xx :



CS-1600/CS-2100/CS-2500 :



CS-5100:

