


# ID-Card Centrifuge

New **HW** ID-Card centrifuge is widely used in hospital, clinical labs and blood donation centers. This is specially centrifuge developed for ID-System with a capacity of 24 ID-Cards.

## PRODUCT FEATURES:

- Maintenance free brushless motor
- LCD display of time and speed
- Fixed speed 870r/min
- Timer from 15s to 99 min. or continuous operation
- Suitable for 24 ID-Cards
- Lid automatically open at the end of the run to prevent samole warming and to reduce handig time



CE compliant according to IVD Directive 98/79/EC

# Centrifuga ID-Card

Nuova Centrifuga per la determinazione su Schedina ID-Card GEL dei gruppi ABO, RhD, CDE e COOMBS in un'unica procedura.

## SPECIFICHE:

- Motore esente da manutenzione
- Display LCD
- Velocità 870 rpm
- Impostazione timer da 15 a 99 min. o funzionamento in continuo
- 24 posizioni per ID-Cards
- Apertura automatica del coperchio a fine lavoro per evitare il surriscaldamento delle ID-Card



### Ordering Information:

**ID-Card Centrifuge**

Product Code: C-17000

Packaging: 1

Meridian Healthcare®

Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 725 6869 (0039)  
Fax: 095 725 4454 (0039)