


## LISS COOMBS Indirect and direct antiglobulin test

ID-Card "LISS/Coombs" with 6 microtubes containing polyspecific AHG (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix.

- Simple, Accurate, Fast.
- Blood samples into citrate, EDTA or CPD-A
- 🕒 Reliable Results in 10 minutes.
- microtubes containing polyspecific AHG
- Direct antiglobulin test (DAT)
- Antibody screening (IAT)
- Compatibility test
- Control



## LISS COOMBS test diretto e indiretto dell'Antiglobulina

Reagenti polispecifici anti-globulina umana (AHG), utilizzati di routine per l'individuazione e l'identificazione di alloanticorpi per test di compatibilità e per il test diretto all'antiglobulina umana (DAT).

- Semplice, Accurato, Veloce.
- Test diretto su citrato, EDTA o CPD-A.
- Risultato in 10 minuti.
- Microprovette contenenti AHG polispecifico in matrice di GEL.
- Test diretto all'antiglobulina (DAT)
- Screening anticorpi (IAT)
- Test di compatibilità
- Controllo



### Ordering Information:

#### Scheda-ID LISS/Coombs

Product Code: 004014V

Packaging: 4x12 pcs



Meridian Healthcare srl  
68 - Tremestieri Etneo  
Catania, IT 95030  
Phone: 095 725 6869 (0039)  
Fax: 095 725 4454 (0039)