

VACUUM TUBES for ESR SEDIMENTATION

Vacuum glass tubes containing Sodium Citrate (3,2% - 0,105M) as anticoagulant



- At least 15 month expiry date
- Aspiration volume 1.28 mL.
- Glass tubes 8×120 mm
- Black butyl rubber cap
- E-Beam Irradiated



Ordering information

REF	PRODUCT	PACK-SIZE
15007	Vacuum Tubes for ESR	100 pcs.

PROVETTE VES per MICROSED-SYSTEM®

Provette in vetro contenente Sodio Citrato (3,2% - 0,105M) come anticoagulante

- Provette in vetro per analizzatori **MICROSED-SYSTEM®**
- Irradiate con tecnologia E-Beam
- Volume di aspirazione: 1.28 mL.
- Tappo in gomma butilica nera
- Provetta in vetro 8×120 mm
- Scadenza minima: 15 mesi

