

Mononucleosis (IM) Rapid Test

The MONO Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of Infectious Mononucleosis heterophile antibodies in whole blood, serum or plasma as an aid in the diagnosis of Infectious Mononucleosis.


- Direct Testing in serum, plasma and wb.
- Fingerstick Whole Blood
- Sample: 25µL. s/p and 50µL. wb
- ⌚ Reliable Results in 5 minutes.
- Sensitivity: 99,9%
- Specificity: 98,9%
- Negative Control: Provided
- Positive Control: Provided



Mononucleosi (IM) Test Rapido

Test rapido per la diagnosi di Infezione da Mononucleosi Infettiva (IM). Rilevazione qualitativa degli anticorpi eterofili anti Mononucleosi Infettiva su sangue intero, siero o plasma. Solo per uso diagnostico professionale in vitro.

- Semplice, Accurato, Veloce.
- Test diretto su sangue intero, siero o plasma.
- Campione: 25µL. per s/p e 50µL per sangue capillare.
- ⌚ Risultato in 5 minuti.
- Sensibilità: 99,9%
- Specificità: 98,9%
- Controlli Positivo e Negativo inclusi.



ORDERING INFORMATION

Codice Dispositivo
IMO-402

Dispositivo
MONONUCLEOSIS (IM) Rapid Test

Confezione
25 test