

# NOROVIRUS Rapid Test

Rapid, qualitative, lateral flow immunoassay test for the detection of Norovirus Antigens Faeces. The assay uses Genogroup 1 and Genogroup 2 specific monoclonal antibodies coated on the test membrane.


- ❖ Simple, Accurate, Fast.
- ❖ Direct Testing on faeces
- ❖ Simultaneous detection of Genogroups I and II
- 🕒 Reliable Results in 15 minutes.
- ❖ Sensitivity : 95,7%
- ❖ Specificity: 91,7%



# NOROVIRUS Rapid Test

Test rapido qualitativo a lettura visiva per la determinazione del NOROVIRUS Antigene su campioni fecali. Il test utilizza gli anticorpi monoclonali specifici di Genogruppo 1 e Genogruppo 2.

- ❖ Semplice, Accurato, Veloce.
- ❖ Test Diretto su campioni fecali.
- ❖ Determinazione simultanea del Genogruppo I e II.
- 🕒 Risultato in 15 minuti.
- ❖ Sensibilità : 95,7%
- ❖ Specificità: 91,7%



## ORDERING INFORMATION:

**Code No.**  
**INO-602**

**Product**  
**NOROVIRUS Rapid Test**

**Packaging**  
**10 Test**