

## NEISSERIA GONORRHOEAE Rapid Test

Neisseria gonorrhoeae Antigen Rapid Test is an immunochromatographic assay for the qualitative presumptive detection of *Neisseria gonorrhoeae* in female endocervical swab and male urethral swab specimens. This kit is intended for use as an aid in the diagnosis of Gonorrhoea.

### FEATURES



- Simple, Accurate, Fast.
- Specimens: endocervical swab and male urethral swab
- 🕒 Reliable Results in 15 minutes.
- Endocervical Sensitivity : 96.9%
- Endocervical Specificity: 95.8%
  - Overall Agreement: 95.9%
- Urethral Sensitivity : 97.8%
- Urethral Specificity: 99.5%
- Overall Agreement: 99.1%

### ORDERING INFORMATION:



Product Code	Product	Pack-Size
IGO-502	Neisseria Gonorrhoeae Antigene Rapid Test	20 test

## NEISSERIA GONORRHOEAE Test Rapido

Neisseria Gonorrhoeae Antigene è un test immunocromatografico per la rilevazione rapida e qualitativa della presenza di *Neisseria gonorrhoeae* da tampone endocervicale femminile e da tampone uretrale maschile. Questo kit è destinato ad essere utilizzato come ausilio nella diagnosi dell'infezione da Gonorrea.

### SPECIFICHE

- Semplice, Accurato, Veloce.
- Campione: tampone endocervicale e uretrale
- 🕒 Risultato in 15 minuti.
- Sensibilità Endocervicale: 96.9%
- Specificità Endocervicale: 95.8%
- Overall Agreement: 95.9%
- Sensibilità Uretrale : 97.8%
- Specificità Uretrale: 99.5%
- Overall Agreement: 99.1%

