


Total Prostate Specific Antigen (t-PSA) Rapid Test

One Step Immunochromatographic assay to detect total prostate specific antigen (PSA) in human serum
For professional use only.


- Simple, Accurate, Fast.
- Sample: human serum.
- ⊕ Reliable Results in 5 minutes.
- Sensitivity : > 98.0%
- Specificity: > 98.0%
- Cut-Off: 4 ng/ml



Total Prostate Specific Antigen (t-PSA) Rapid Test

Test rapido immunocromatografico Lateral Flow per la rilevazione dell'antigene specifico prostatico in campioni di siero umano.

- Semplice, Accurato, Veloce.
- Test diretto su siero umano.
- ⊕ Risultato in 5 minuti.
- Sensibilità : > 98.0%
- Specificità : > 98.0%
- Cut-Off: 4 ng/ml.



EUROPE Ordering Information:

PSA RAPID TEST
Product Code: 1504021
Packaging: 20 test

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)