

## NORO+ROTA+ADENOVIRUS Rapid Test

The Norovirus/Rotavirus/Adenovirus Combo Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of norovirus, rotavirus and adenovirus in human feces specimens to aid in the diagnosis of norovirus, rotavirus or adenovirus infection.


- ❖ Simple, Accurate, Fast.
- ❖ Direct Testing on faeces
- ❖ Simultaneous detection of NORO+ROTA+ADENO
- 🕒 Reliable Results in 15 minutes.
- ❖ NORO Sensitivity : 95,7%
- ❖ NORO Specificity: 91,7%
- ❖ ROTA Sensitivity: 97,3%
- ❖ ROTA Specificity: 97,1%
- ❖ ADENO Sensitivity: 95,2%
- ❖ ADENO Specificity: 97,7%



## NORO+ROTA+ADENOVIRUS Test Rapido

Test rapido qualitativo a lettura visiva per la rilevazione differenziata di ADENOVIRUS/ROTAVIRUS/NOROVIRUS Antigeni su campioni fecali. Il NORO utilizza anticorpi monoclonali specifici di Genogruppo 1 e Genogruppo 2.

- ❖ Semplice, Accurato, Veloce.
- ❖ Test Diretto su campioni fecali.
- ❖ Rilevazione simultanea di NORO+ROTA+ADENO.
- 🕒 Risultato in 15 minuti.
- ❖ NORO Sensibilità : 95,7%
- ❖ NORO Specificità: 91,7%
- ❖ ROTA Sensibilità: 97,3%
- ❖ ROTA Specificità: 97,1%
- ❖ ADENO Sensibilità: 95,2%
- ❖ ADENO Specificità: 97,7%



### ORDERING INFORMATION:

**Code No.**  
**INRA-635**

**Product**  
**NORO+ROTA+ADENO Rapid Test**

**Packaging**  
**25 Test**