

MANUFACTURE STANDARD, MANUFACTURE RELIANCE

RT-A29

Microplate Incubator & Shaker

COMPACT SIZE
USE FOR ELISA READER
OR CELL CULTURE PLATE



Features

- Compatible for Elisa Microplate or Cell Culture Plate
- Incubates range from 25°C to 60°C
- Monochrome screen display
- 250-1200 rotation per minute for frequency



RT-A29

Microplate Incubator & Shaker



Technical Specifications

Power Supply	AC 100V-240V, 50Hz/60Hz
Temperature Range	25-60°C
Temperature Accuracy	±0.5°C
Temperature Uniformity	250-1200rpm (RPM= rotation per minute)
Amplitude	2mm
Timin	0-100 hours
Compability	Compatible for Elisa Microplate or Cell Culture Plate
Temperature Display	Monochrome screen
Dimension L×W×H(mm)	270mm*270mm*130mm

Meridian Healthcare srl

Tel. +39 095 725 68 69 Fax: +39 095 725 44 54
info@meridianhealthcare.it
www.meridianhealthcare.it

Meridian Healthcare®



Meridian Healthcare®

Trade diagnostical products
for clinical laboratories
products for life.

