

Rickettsia IgG IgM Combo Rapid Test

Rapid Rickettsia IgG/IgM Combo Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Rickettsia in human whole blood, serum or plasma. The assay is used as a screening test for Rickettsia infection and as an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria. The Rickettsia test is a laboratory test used to detect human antibodies to Rickettsia bacteria in a blood sample. Rickettsia rickettsii is the causative agent of Rocky Mountain spotted fever. The illness may vary greatly in severity and presentation from very mild illness through to life-threatening cases. If diagnosed early rickettsial diseases are readily treatable with antibiotics

IgG AND IgM ANTIBODIES TO RICKETTSIA

FEATURES

- Simple, Accurate, Fast.
- Sample: serum, plasma, wb.
- 🕒 Reliable Results in 20 minutes.
- Combo IgG/IgM
- Sensitivity: 100%
- Specificity: 100%
- CE marked



RMSF infectious

Rickettsia IgG IgM Combo Rapid Test

Il test rapido combo IgG / IgM Rickettsia è un dispositivo immunocromatografico per la rilevazione simultanea di anticorpi IgG e IgM contro Rickettsia su sangue intero, siero o plasma umani. Il test viene utilizzato come screening per l'infezione da Rickettsia e come ausilio per la diagnosi differenziale delle infezioni primarie e secondarie in combinazione con altri criteri.

Il test Rickettsia è un test di laboratorio utilizzato per rilevare anticorpi umani contro i batteri Rickettsia su campioni di sangue. Rickettsia rickettsii è l'agente causale della febbre maculosa di Rocky Mountain. La malattia può variare notevolmente in gravità e presentazione da malattia molto lieve a casi potenzialmente letali. Se diagnosticato rapidamente, le infezioni da rickettsie sono facilmente curabili con antibiotici.

SPECIFICHE

- Semplice, Accurato, Veloce.
- Campione: siero, plasma e sangue intero.
- Combo IgG e IgM
- 🕒 Risultato in 20 minuti.
- Sensibilità: 100%
- Specificità: 100%
- CE marked



EUROPE Ordering Information:



PRODUCT

Rickettsia IgG/IgM Combo (71017 U.S.A.)

REF.

1006C2

N. Test

25 card