

STREPTOCOCCUS β -hemolytic group A Rapid Test

Strep A test is a lateral flow, one-step immunoassay for the rapid qualitative detection of Group A Streptococcal antigen directly from throat swabs. Only for laboratory use.



- Simple, Accurate, Fast
- Direct Testing **throat swabs**
- Ambulatory can be used
- 🕒 Reliable Results in 5 minutes
 - Sensitivity : >97%
 - Specificity: >97%
- Positive Control: 1 vial
 - CE marked

STREPTOCOCCO β -emolitico gruppo A Rapid Test

Test rapido immunocromatografico Lateral Flow per la determinazione qualitativa dell'antigene dello streptococco di gruppo A in tamponi faringei.

- Semplice, Accurato, Veloce
- Test diretto da tampone faringeo
- Utilizzo in laboratorio e in ambulatorio
- 🕒 Risultato in 5 minuti
- Sensibilità : >97%
- Specificità: >97%
- Controlli inclusi nel kit
- CE marked

Controlli terzi disponibili:

- Ctrl Positivo: ODZ-212 kit da 10 test
- Ctrl Positivo: RM07110PC fiala da 1 ml



EUROPE Ordering Information:

STREPTO β -hemolytic group A RAPID TEST

Product Code: 222001A

Packaging: 20 test + ctrl POS.



Product Code: ISTB-502

Packaging: 20 test + ctrl POS. - NEG.

Meridian Healthcare®

Meridian Healthcare srl
68 - Tremestieri Etneo
Catania, IT 95030
Phone: 095 725 6869 (0039)
Fax: 095 725 4454 (0039)