

SYPHILIS Treponema Pallidum Rapid Test

The SYPHILIS Test Device is a rapid chromatographic immunoassay for the qualitative detection of antibodies IgG and IgM to Treponema Pallidum in serum, plasma or whole blood.



- Simple, Accurate, Fast.
- Direct Testing on serum, plasma and whole blood.
- 🕒 Reliable Results in 10 minutes.
- Sensitivity : > 99.9%
- Specificity: 99.7%

ORDERING INFORMATION:



Product Code
ISY-402

Product
SYPHILIS Ab Rapid Test

Pack-Size
25 test

SYPHILIS Treponema Pallidum Rapid Test

Test rapido immunocromatografico Lateral Flow per la rilevazione qualitativa degli anticorpi IgG e IgM anti-Treponema Pallidum su campioni di siero, plasma e sangue intero.

- Semplice, Accurato, Veloce.
- Test diretto su siero, plasma e sangue intero
- 🕒 Risultato in 10 minuti.
- Sensibilità : > 99.9%
- Specificità: 99.7%

JusChek

