


## SCANGEL LINE ABO - RH ID-Card DiaClon

The ID-Card "DiaClon ABO/Rh for Patients" offers a complete test profile for forward typing of ABO plus RhD and CDE status in one single test procedure.

- Simple, Accurate, Fast.
- Blood samples into citrate, EDTA or CPD-A
- ⌚ Reliable Results in 10 minutes.

ID-Card "DiaClon ABO/Rh for Patients" contains:

- monoclonal anti-A [cell line A5]
- anti-B [cell line G½]
- anti-AB [cell lines ES 131 (ES-15), Birma-1, ES-4]
- anti-D [cell lines LHM 59/20 (LDM 3) and 175-2]
- anti-CDE (cell lines MS-24, MS-201, MS-26, MS-80)
- Control




## SCANGEL LINE ABO - RH scheda-ID DiaClon

La scheda-ID "DiaClon ABO/Rh for Patients" consente di ottenere un profilo di analisi completo per ABO, RhD e CDE in un'unica procedura.

- Semplice, Accurato, Veloce.
- Test diretto su emazie diluite.
- ⌚ Risultato in 10 minuti.

La scheda ID DiaClon contiene nella matrice di Gel i seguenti antisieri:

- anti-A monoclonale [linea cellulare: A5]
- anti-B [linea cellulare: G½]
- anti-AB [linee cellulari: ES 131 (ES-15), Birma-1, ES-4]
- anti-D [linee cellulari: LHM 59/20 (LDM 3) e 175-2]
- anti-CDE (linee cellulari: MS-24, MS-201, MS-26, MS-80)
- Controllo



Ordering Information:

**DiaClon ABO/Rh**

Product Code: 001044V

Packaging: 4x12 pcs ( 48 test )

Meridian Healthcare®

Meridian Healthcare srl

68 - Tremestieri Etneo

Catania, IT 95030

Phone: 095 725 6869 (0039)

Fax: 095 725 4454 (0039)