



# TRICHOMONAS Vaginalis Rapid Test

New Rapid Test Lateral Flow for the qualitative detection of TRICHOMONAS Vaginalis n swab from vaginal cavity and URINA.

- Simple, Accurate, Fast
- Sample: vaginal swab and URINA
- ⌚ Reliable Results in 15 minutes
- Sensitivity: 100%
- Specificity: 99%

### PERFORMANCE

The gold standard was the patient positive to speculum and vaginal examinations to be diagnosed Trichomonas Vaginalis infection. The total 204 cases were compared to wet mount microscopy pouch culture and JD's TVAG Kit. The JD's TVAG test kit will detect soluble antigen present in vaginal samples from > 20 organisms, a number of organisms that is lower than that expected from discharge of patients.



### ORDERING INFORMATION:

**Product Code**  
JTVAG-02

**Product**  
TRICHOMONAS Vaginalis Rapid Test

**Pack Size**  
25 test

