

# anti-CCP

Anti-Cyclic Citrullinated Peptide  
semi-quantitative immunochromatographic test

## LATERAL FLOW

Semi-quantitative Anti-CCP measurement  
Rapid, accurate, lab quality results within minutes  
All-inclusive kit including reference card for value estimation  
Limit of detection: 25 RU/ml.  
Sensitivity: 90%  
Specificity: 99,4%



# Anti-cyclic Citrullinated Peptide (anti-CCP) semi-quantitative assay kit Colloidal Gold

This kit is used for semi-quantitative measurement of human anti-cyclic citrullinated peptide antibody (Anti-CCP) in serum and plasma employing colloidal gold immunochromatography as an aiding tool for diagnosis of rheumatoid arthritis (RA).

**Negative**      **Positive**      **Invalid**

Anti-Cyclic Citrullinated Peptide Antibody (Anti-CCP)      Test Result (RU/mL)

Patient Name \_\_\_\_\_

Sample ID \_\_\_\_\_

Operator \_\_\_\_\_

Test Time \_\_\_\_\_

Control Line (C) →

Test Line (T) →

<input type="checkbox"/>	>400	
<input type="checkbox"/>	400	400
<input type="checkbox"/>	250	
<input type="checkbox"/>	100	100
<input type="checkbox"/>	62.5	
<input type="checkbox"/>	25	25
<input type="checkbox"/>	<25	

**ORDERING INFORMATION:**

<b>Product Code</b> HGCG02	<b>Items</b> Anti-CCP Test Rapido	<b>Packaging</b> 40 test
-------------------------------	--------------------------------------	-----------------------------