

# TUBERCULOSIS IgG IgM Combo Rapid Test

The TB IgG IgM Rapid Test is a sandwich lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of IgM anti-Myco bacterium Tuberculosis and IgG anti-M.TB in human whole blood, serum or plasma.

## FEATURES

- Simple, Accurate, Fast.
- Direct Testing on serum, plasma or whole blood.
- ⌚ Reliable Results in 15 minutes
- Sensitivity IgG : 88,6%
- Specificity IgG: 96,5%
- Sensitivity IgM : 85,7%
- Specificity IgM: 96,5%

JusChek®



# TUBERCOLOSI IgG IgM Combo Test Rapido

Test rapido immunocromatografico lateral flow per la rilevazione qualitativa e differenziata degli anticorpi anti-M TB IgG e IgM su campioni di siero, plasma e sangue intero.



## SPECIFICHE

- Semplice, Accurato, Veloce.
- Test diretto su sangue intero.
- Utilizzo previsto anche in ambulatorio.
- ⌚ Risultato in 15 minuti.
- Sensibilità IgG: 88,6%
- Specificità IgG: 96,5%
- Sensibilità IgM: 87,5%
- Specificità IgM: 96,5%
- CE marked

Negative

Invalid

Positive



IgG  
Positive

IgM  
Positive

IgG/IgM  
Positive

EUROPE Ordering Information:



### PRODUCT

TUBERCULOSIS IgG IgM Combo

### REF.

ITB-T402

### N. Test

25 card