


BIOPANDA™ MEASLES IgG/IgM Rapid Test

The Biopanda Measles IgG/IgM Rapid Test qualitatively detects IgG and IgM antibodies to measles in human whole blood, serum, or plasma samples. This test applies lateral flow immuno-chromatography and aids in the diagnosis of measles infections.



- Direct Testing in serum, plasma and wb.
- Sample: 10µL. s/p and 20µL. wb
- 🕒 **Reliable Results in 15 minutes.**
- Sensitivity IgG : 97,9%
- Specificity IgG: 95,7%
- Sensitivity IgM : 94,7%
- Specificity IgM: 94,1%



BIOPANDA™ MORBILLO IgG/IgM Test Rapido

Il dispositivo **BIOPANDA™ MEASLES** è un test rapido Lateral Flow per la determinazione simultanea e differenziata degli anticorpi anti-Morbillo IgG e IgM su campioni di sangue intero, siero o plasma.

- Semplice, Accurato, Veloce.
- Test diretto su sangue intero, siero o plasma
- Campione: 10µL. per s/p e 20µL per sangue venoso o capillare.
- 🕒 **Risultato in 15 minuti.**
- Sensibilità IgG: 97,9%
- Specificità IgG: 95,7%
- Sensibilità IgM: 94,7%
- Specificità IgM: 94,1%



ORDERING INFORMATION

Codice Dispositivo
RAPG-MSL-001

Dispositivo
BIOPANDA™ Measles IgG/IgM Rapid Test

Confezione
20 test