

C1 ESTERASE INHIBITOR TURBI IMMUNOASSAY

Determination of the C1 esterase inhibitor based on an immuno-turbidimetric antigen-antibody reaction.
For *in vitro* diagnostic use only.



- Immuno-turbidimetric antigen-antibody reaction.
 - Dual Vial Liquid Stable Test Kit .
 - 🕒 Reliable Results in 5 minutes.
 - Measuring range: 0 - 90 mg/dl.
 - Sample: Serum.
- Applications are available for many analyzers after request.

For calibration use
Centronic Protein Calibration Set (PROCAL).



C1 ESTERASI INIBITORE TURBI IMMUNOASSAY

Saggio immuno-turbidimetrico per la determinazione della proteina C1 Esterasi Inibitore (C1-INH) su campioni di siero.

- Metodo Immuno-turbidimetrico.
- Due reagenti liquidi.
- 🕒 Risultato in 5 minuti.
- Saggio su siero.
- Automatizzabile su tutte le piattaforme open system.
- Range di misura: 0 - 90 mg/dl.



ORDERING INFORMATION

REF.

CF19000055
PROCAL

PRODUCTS

C1 Esterase Inhibitor
Protein Calibration Set

PACKAGING

2x25ml/1x5ml
5x1 ml