

# ALLTEST® HIV 1/2 Rapid Test

The ALLTEST® HIV Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies to Human Immunodeficiency Virus (HIV) type 1 (including subtype O) and type 2 in whole blood, serum or plasma to aid in the diagnosis of an HIV infection.

- ✓ Simple, Accurate, Fast.
- ✓ Direct Testing in serum, plasma or whole blood.
- 🕒 Reliable Results in 10 minutes.
- ✓ HIV-1 sensitivity : >99%
- ✓ HIV-2 sensitivity : >99%
- ✓ Specificity: 100%



CE 2934

## ORDERING INFORMATION:

Code No.	Product	Pack-size
IHI-402	HIV ½ Rapid Test	25 Test

# ALLTEST® HIV 1/2 Test Rapido

Test rapido per la rilevazione degli anticorpi contro l'HIV (virus dell'immunodeficienza umana) di tipo 1 (compreso il sottotipo O) e di tipo 2 su campioni di sangue intero, siero o plasma.

- ✓ Semplice, Accurato, Veloce.
- ✓ Test diretto su siero, plasma o sangue intero.
- 🕒 Risultato in 10 minuti.
- ✓ Sensibilità: HIV-1 >99%.
- ✓ Sensibilità: HIV-2 >99%.
- ✓ Specificità: 100%.

